

MadaCide-FDW-PLUS

No. 7032 Disinfecting/Deodorizing/Cleaning Surface Wipes

EFFECTIVE AGAINST:

- ✓ Staphylococcus aureus
- ✓ Pseudomonas aeruginosa
- ✓ Salmonella enterica
- ✓ Mycobacterium tuberculosis
- ✓ Methicillin Resistant Staphylococcus aureus - MRSA
- ✓ Vancomycin Resistant Enterococcus faecalis - VRE
- ✓ Influenza A - (A2 Japan)
- ✓ Escherichia coli O157:H7
- ✓ Human Hepatitis B virus
- ✓ Human Hepatitis C virus
- ✓ Herpes Simplex virus 2 (genital herpes virus)
- ✓ HIV-1 (AIDS virus) (At an exposure time of 1 min)
- ✓ Trichophyton mentagrophytes (Athlete's foot fungus)



MadaCide-FDW-PLUS DISINFECTING/DEODORIZING/CLEANING WIPES

FOR USE ON ENVIRONMENTAL SURFACES

ACTIVE INGREDIENTS:

Virucidal*	N-Alkyl (60% C14, 30% C16, 5% C12, 5% C18).....0.12%
	Dimethyl Benzyl Ammonium Chloride
Bactericidal	N-Alkyl (68% C12, 32% C14).....0.12%
	Dimethyl Ethyl Benzyl Ammonium Chloride
Pseudomonocidal	Isopropyl Alcohol.....41.58%
Tuberculocidal	
Fungicidal	
	OTHER INGREDIENTS:.....58.18%
	TOTAL:.....100.00%

* Does not include weight of dry wipe

KEEP OUT OF REACH OF CHILDREN WARNING

DO NOT USE AS A BABY WIPE

See back panel for additional precautionary statements.
A non-woven disposable cloth containing a stable, low pH formulated disinfectant and deodorant for use on hard, non-porous surfaces where control of cross-contamination is required.

Reorder No. 7032

NET CONTENTS: 160 Pre-saturated Wipes (6x6.75 in)
NET WEIGHT: 1lb. 10 ounces

MadaCide-FDW-PLUS Disinfecting/Deodorizing/Cleaning Wipes are effective against the following pathogens in 2 min at 20°C/68°F, on hard, non-porous surfaces when used:

Staphylococcus aureus Pseudomonas aeruginosa Salmonella enterica
Mycobacterium tuberculosis Methicillin Resistant Staphylococcus aureus - MRSA
Escherichia coli O157:H7 Vancomycin Resistant Enterococcus faecalis - VRE
Human Hepatitis B virus Human Hepatitis C virus HIV-1 (AIDS virus)
*See back panel for complete list of pathogens

Kills Pandemic 2009 H1N1 Influenza A virus (formerly swine flu)

MadaCide-FDW-PLUS Disinfecting/Deodorizing/Cleaning Wipes are effective against the following pathogens in 5 minutes at 20°C/68°F, on hard, non-porous surfaces when used as directed:
Trichophyton mentagrophytes (Athlete's foot fungus)

Areas of Use:
hospital emergency medical units
intensive care units emergency rooms
ambulances emergency medical services
nursing homes nursing facilities
orthopedics respiratory therapy
patient rooms health care
medical offices physician offices
dental offices exam rooms
veterinary facilities lockers
janitoria other similar areas

PREPARATION OF DISPENSER: Remove cover and discard covering seal. From center of wiper roll, pull wiper through the lid cover, twist into a panel and thread through the top hook. Pull wiper through about one (1) inch on container, ensuring a secure fit. Pull out lid and seal on container, ensuring a secure fit. All out lid and seal at angle. Remaining wipes feed automatically for next use. Hinged cap after each use to insure moisture retention.

160 Pre-Saturated Wipes (6" x 6.75") 12/case

EPA REGISTERED



625 WASHINGTON AVE, CARLSTADT, NJ 07072
CALL 1-800-526-6370 | FAX 201-460-3509 | WEB www.madamedical.com

MadaCide-FDW-PLUS

MadaCide-FDW-PLUS Wipes are non-woven disposable cloths containing a stable, low pH formulated disinfectant and deodorant cleaner. They are designed for use on hard, non-porous surfaces.

ORGANISM(S)

CONTACT TIME

Staphylococcus aureus, Pseudomonas aeruginosa, Salmonella enterica, Mycobacterium tuberculosis, Methicillin Resistant Staphylococcus aureus - MRSA, Vancomycin Resistant Enterococcus faecalis - VRE, Influenza A - (A2 Japan), Escherichia coli 0157:H7, Human Hepatitis B virus, Human Hepatitis C virus, Herpes Simplex virus 2 (genital herpes virus), HIV-1 (AIDS virus)** (At an exposure time of 1 min)

2 min

Trichophyton mentagrophytes (Athlete's foot fungus)

5 min

AREAS OF USE:

Stainless steel, formica, tables, carts, baskets, counters, cabinets, and telephones. This product can be used in hospitals, intensive care units, surgery, recovery, anesthesia, x-ray, catheterization laboratory, orthopedics, newborn nursery, respiratory therapy, emergency medical settings, clinics, nursing facilities, medical offices, dental offices, veterinary, and janitorial facilities where control of cross-contamination is required.