

VivaLIFT!®

POWER RECLINER COLLECTION

THE MOST
COMFORTABLE
RECLINER
IN THE WORLD!

THE 4-MOTOR RECLINER



WITH
COUNTLESS
POSITION
COMBINATIONS

LEGACY

Embrace the serenity of a quiet corner in your home with a cozy, luxurious recliner. Take time for yourself with the **Legacy Collection** from VivaLift!®. Put your feet up and fully stretch out on the footrest extension. Enjoy a good novel or turn on the television. No matter what activity you choose, infinite, lay-flat positioning ensures you'll always find your most comfortable and supportive position. Plus, you'll never misplace your book or the remote with the convenient side storage pocket.









Scan here for
more information
on **VivaLift!®**




PRIDE®
MOBILITY
LIVE YOUR BEST®

COLLECTION	LEGACY	LEGACY
MODEL	PLR-958M	PLR-958L
POSITION	Infinite Lay Flat Relaxed Reclined Upward Lift	Infinite Lay Flat Relaxed Reclined Upward Lift
SIZE	Medium	Large
WEIGHT CAPACITY	400 lbs.	400 lbs.
SEAT-TO-FLOOR	19.5"	20.5"
SEAT DEPTH	21.5"	22"
SEAT WIDTH	22"	23"
TOP OF BACK TO SEAT	25"	26"
FULLY PADDED CHAISE	Yes	Yes
FOOTREST EXTENSION	Standard	Standard
SEATING TYPE	Total Comfort	Total Comfort
CHAIR HEIGHT	42"	43.5"
DISTANCE FROM WALL	16"	16"
OVERALL WIDTH	35"	36.5"
CHAIR WEIGHT	153 lbs.	160 lbs.
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated

INDEPENDENT FUNCTIONS

-  1. The power headrest lifts your head, neck and shoulders for maximum support.
-  2. The upward lift offers convenient stand assist for getting in and out of the recliner.
-  3. The power lumbar enables you to personalize your comfort.
-  4. The footrest extension allows you to relax and stretch out completely.
-  5. The backrest moves up and down independently to help you find the perfect position.
-  6. Infinite position recline enables the backrest and footrest to move independently from one another so that a variety of positions can be reached.



STANDARD FEATURES

7. The USB remote adjusts all functions independently and allows you to charge your smart device.
8. Lithium battery backup ensures your recliner is operational when the power goes out.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

 = Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

STANDARD FABRIC: SAVILLE



WINE

GREY

BROWN

(US) 800-800-8586 / (CAN) 888-570-1113
prideliftchairs.com / pridemobility.com

PRIDE
MOBILITY
LIVE YOUR BEST®