

# VivaLIFT!®

POWER RECLINER COLLECTION

THE MOST  
COMFORTABLE  
RECLINER  
IN THE WORLD!

## THE 4-MOTOR RECLINER



Shown in Astro Brown Fabric

WITH  
COUNTLESS  
POSITION  
COMBINATIONS

## TRANQUIL

There is no greater joy than reading a good book in a comfortable, cozy recliner. The **Tranquil Collection** by VivaLift!® Power Recliners combines quiet comfort with gorgeous fabrics, for a rich and upscale look. Get the support you need for your neck, shoulders and lower back as you dive into your favorite story. The Tranquil Collection exudes style and class with contrast stitching. It's the perfect addition to your home.



Scan here for  
more information  
on **VivaLift!®**









**PRIDE**®  
MOBILITY  
LIVE YOUR BEST®


# TRANQUIL

**VivaLIFT!**  
POWER RECLINER COLLECTION

COLLECTION	TRANQUIL	TRANQUIL	TRANQUIL	TRANQUIL
MODEL	PLR-935PW	PLR-935S	PLR-935M	PLR-935LT
POSITION	Infinite Lay Flat Relaxed Reclined Upward Lift	Infinite Lay Flat Relaxed Reclined Upward Lift	Infinite Lay Flat Relaxed Reclined Upward Lift	Infinite Lay Flat Relaxed Reclined Upward Lift
SIZE	Petite Wide	Small	Medium	Large/Tall
WEIGHT CAPACITY	400 lbs.	400 lbs.	400 lbs.	400 lbs.
SEAT-TO-FLOOR	18"	18"	19.5"	21.5"
SEAT DEPTH	18.5"	18.5"	20.5"	22.5"
SEAT WIDTH	24"	20"	22"	23.5"
TOP OF BACK TO SEAT	24.5"	25"	27"	30"
FULLY PADDED CHAISE	Yes	Yes	Yes	Yes
FOOTREST EXTENSION	Standard	Standard	Standard	Standard
SEATING TYPE	Total Comfort	Total Comfort	Total Comfort	Total Comfort
CHAIR HEIGHT	38"	40"	42"	47"
DISTANCE FROM WALL	12.5"	13"	16"	16"
OVERALL WIDTH	36"	33.75"	37.5"	38"
CHAIR WEIGHT	143 lbs.	137 lbs.	145 lbs.	155 lbs.
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated

## INDEPENDENT FUNCTIONS

-  1. The power headrest lifts your head, neck and shoulders for maximum support.
-  2. The upward lift offers convenient stand assist for getting in and out of the recliner.
-  3. The power lumbar enables you to personalize your comfort.
-  4. The footrest extension allows you to relax and stretch out completely.
-  5. The backrest moves up and down independently to help you find the perfect position.
-  6. Infinite position recline enables the backrest and footrest to move independently from one another so that a variety of positions can be reached.



## STANDARD FEATURES

7. The USB remote adjusts all functions independently and allows you to charge your smart device.
8. Lithium battery backup ensures your recliner is operational when the power goes out.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

 = Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

## STANDARD FABRIC: ASTRO



GREY

BROWN

MUSHROOM

(US) 800-800-8586 / (CAN) 888-570-1113  
prideliftchairs.com / pridemobility.com

**PRIDE**<sup>®</sup>  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>